

Nora Research

Innovative Technology Products

CATALOGUE

05/2012



Nora Research

Professional device for bio-resonance therapy
(Rife/Clark frequency generator)



Heiron



Bio-resonance therapy - alternative

Purpose:

Professional electronic apparatus "Heiron" is a new generation device and is designed for active bio-resonance therapy (ABT) on the human body based on the generation of electromagnetic signals of certain frequencies and waveforms, through which the influence of harmful pathogens (bacteria and viruses) and parasites in the human body and clearing it of toxins. Treatment without medication for a very wide range of diseases.

The device is specifically designed to meet all requirements professional contact bioresonance therapy according to the writings of scientists, Dr. Hulda Clark and Dr. Raymond Rife, and for convenience is attached "Frequency reference guide book" - therapeutic frequencies Rife / Clark. Attached is an additional computer CD-ROM drive with the entire database of frequencies CAFL for over 1,300 diseases. Frequency range is from 0 to 10 000 000 Hz (0-10 MHz) in increments of 1Hz and 0.1Hz accuracy. There are two interchangeable modes according to the type of therapy: a rectangular output signal (for the Clark therapy) and sinusoidal output signal (Rife

therapy). The output power is smoothly adjustable and can be precisely adapted to the sensitivity of each organism and the disease (0-22V). Setting the switching frequencies and modes as easily and are made by pressing quasi-sensor buttons. The output frequency is displayed on the precision illuminated LCD display, which facilitates control in. Available indicator light sensitivity (output control), which can be monitored 100% contact with the body. Built are sound and light control indicators. There are 4 current protection, which ensures complete safety in use. Both types of electrodes (bracelets and handles) kit provides flexibility therapy. We use the most modern technology combined with unique microprocessor program allows the device to have high technical parameters combined with maximum easier use. These properties allow it to be applied by all medical professionals and therapists in their practice, but also for private everyday use. Belongs to a professional grade appliances. Designed and manufactured entirely in the European Union. Available detailed instructions in the English language. There is warranty and post warranty maintenance. More information available on our website:


www.noraresearch.com.

Technical characteristics:

- Spectrum:
 - Mode "Rife": sinusoidal output signal (THD: <0.01%)
 - Mode "Clark": rectangular output signal with duty cycle 50%
- Frequency range:
 - Mode "Rife": 0-9999999 Hz (sine wave)
 - Mode "Clark": 0-7000000 Hz (rectangle)
- Accuracy of the generated frequencies: <0.1Hz
- Resolution of frequency: 1Hz
- Visualization of frequency: STN LED LCD display
- Output Offset: 100% positive
- Output: smoothly adjustable 0 - 22V + / - 10% (5-4000000 Hz)
- Accuracy of marking on the output regulator: <20% (5-4000000 Hz)
- Maximum output current: 0.01A
- Protections: 3 built-in current protections + external removable fuse.
- Power supply unit: a network 220/230V / 50Hz (+ / - 15%)
- Power consumption: max.10W
- Temperature range for normal operation: from +10 to +35 degrees Celsius
- Operation in humidity: <70%
- Weight: not more than 3 kg.

Features

Modern design, comfortable control and display



Reliable and comfortable electrodes



Rich accessories, detailed documentation



www.noraresearch.com
nora.research@mail.bg



Nora Research

Professional quantum therapy device



Hera



Quanta - the energy in the universe

Purpose:

Professional electronic therapeutic device "Hera" is a new generation instrument is designed to perform quantum therapy on the human organism. Quantum Therapy is a modern trend in modern medicine. Used in more than 9500 hospitals in Russia, Europe and USA. Has been shown to decrease the sick patients from 3 to 6 times in Quantum treatment therapy. Impact on biologically active zones and bodies in the human body via electromagnetic radiation within a certain range. Thereby stimulating its own vitality to fight many diseases, improves circulation, activates the body's immune system and has anti-inflammatory and analgesic action. With his help to carry out without drug and painless treatment for a wide range of diseases. Upon acceptance of medication therapy apparatus fold increase in their effectiveness. Unlike medicines, apparatus "Hera" not just treat the disease, but weakened and transformed cancer cells of the organism in sound health, thereby

achieving a strong effect of treatment for a short period of time. Thanks to new technologies and the element base has been possible to have a handset offering therapy, which until recently offered only clinics. There are 4 programs to therapy, and 3 automatic interval duration of therapy. Placing the device and work are most facilitated and carried out by pressing buttons. In combination with modern design provided sound and light indicators. Device is characterized by its comfort, quality and long life, which makes it suitable for use in everyday life, medical and therapeutic centers. Therapeutic head is compact and flexible cable which facilitates work with it. Power is carried through safe transformer adapter. Modern microelectronics is used as the apparatus is designed and manufactured entirely in the European Union. Available is a detailed instruction manual in English. There is a guarantee and warranty support. More information can be found on our website: www.noraresearch.com

Technical characteristics:

Types of radiation:

- Pulse infrared laser spectrum close
- Pulse broad infrared radiation
- Pulsating radiation in the red (visible) spectrum
- Permanently radiation in red (visible) range

Length of rays:

- Laser radiation: 800 to 910 nm
- Infrared radiation: 860 to 960 nm
- Red radiation: 500 to 700 nm

Pulse power of laser infrared radiation: max. 8W

Intensity of radiation in the red spectrum: max. 27,000 mcd

Power transmission in the red spectrum: max. 1W

Time the session: automatic - 5, 10 or 20 min + / - 10%

Power of the device: through transformer power adapter from the network 220/230V / 50Hz (optional 110V)

Power Consumption: max.10W

Mass: no more than 2 kg.


Warranty: 12 months

Features

convenient therapeutic head



4 programs, 3 time intervals, suitable indicators and buttons



modern and comfortable design



www.noraresearch.com
nora.research@mail.bg



Nora Research

Professional electronic device for bio-resonance therapy
(professional Zapper)



Tangra II



Bio-resonance therapy - alternative

Purpose:

Professional electronic device "Tangra 2" is a new generation instrument is designed for active bio-resonance therapy (ABT) on the human body based on the generation of electromagnetic signals with certain frequencies and form by which the influence of harmful pathogens (viruses and bacteria) and parasites in human body. Typical for this type of therapy that is totally harmless, and that helps in the case where no other method does not help. Daily treatment with device clears the body of harmful parasites that can cause hundreds of diseases, while increasing tone, removes fatigue and improves overall health.

It is based entirely on the writings of the famous Dr. Hulda Clark, which laid the foundations of modern bioresonance therapy practice and through your great books ("Treatment of all diseases, etc..). There are 4 work therapy programs, which allows to have an influence on a wide range of harmful

microorganisms. One of the programs is that of electronic Zapper Dr. Hulda Clark. Available is a healing program for high-frequency therapy, which acts as an immune stimulant. Possible change in the power of the device, thereby achieving a better setting to the disease and given the sensitivity of human organism.

Device is characterized by its comfort, quality and long life, which makes it suitable for use in everyday life, medical and therapeutic centers. There are healthy and convenient removable bracelets electrodes. Power is carried through safe transformer adapter. Contemporary design, comfortable buttons, light and sound indicators make use of the device comfortably. Modern microelectronics is used as the unit. Is designed and manufactured entirely in the European Union. Available is a detailed instruction manual in English. There is a guarantee and warranty support. More information can be found on our website: www.noraresearch.com

Technical characteristics:

Power supply: 220V/50Hz mains via transformer adapter

Power consumption: <2W

Spectrum:

Program I: Time and frequency modulated signal with a bearing rate 33 000 Hz - rectangular shape with a wide harmonic range

Program II: Time and frequency modulated signal with a bearing rate 1000 Hz - rectangular shape with a wide harmonic range

Program III: Time and frequency modulated signal with a bearing rate 5000 Hz - rectangular shape with a wide harmonic range

Program Ultra: Time and frequency modulated signal with high frequency (HFT therapy) - a rectangular shape with a wide harmonic range

Amplitude of output signal: 25-30V

Maximum output current:

Mode "Low" 80 micro amps


Mode "High" 200 micro amps

Time at work: 7 minutes + / -30 seconds

Warranty: 12 months

Features

convenient removable bracelets electrodes



4 programs, 2 power levels
suitable indicators and buttons



modern and comfortable design



www.noraresearch.com
nora.research@mail.bg



Nora Research

Electronic device for bio-resonance therapy
(Zapper)



Tangra



Bio-resonance therapy - alternative

Purpose:

Professional electronic device "Tangra" is a new generation instrument is designed for an active bio-resonance therapy (ABT) on the human body based on the generation of electromagnetic signals with certain frequencies and form through which the harmful influence pathogenic (bacteria and viruses) and parasites in the human body. Typical for this type of therapy that is totally harmless, and that helps in the case where no other method does not help. Daily treatment with apparatus clears the body of harmful parasites that can cause hundreds of diseases, while increasing tone, removes fatigue and improves overall health. It is based entirely on the writings of the famous Dr. Hulda Clark, which laid the

foundations of modern bio-resonance therapy practice and through his great books ("Treatment of all diseases, etc.). Device "Tangra" is characterized by its comfort, quality and long life, which makes it suitable for use in everyday life, medical and therapeutic centers. There are comfortable and healthy bracelets electrodes that make carrying comfortable therapy. Device can be powered by battery or by network, making it independent of the situation. Modern microelectronics is used as the device is designed and manufactured entirely in the European Union. Available is a detailed instruction manual in English. There is a guarantee and warranty support. More information can be found on our website:

www.noraresearch.com

Technical characteristics:

Power Supply: Mains 220V/50Hz (110V U.S. version) by transformer adapter or 9V battery type 6F22

Power consumption: <0.5W

Spectrum: Time and frequency modulated rectangular signal with a wide harmonic spectrum consistent with the scientific writings of Dr. Hulda Clark.

Amplitude of the output: max. 30V


Maximum output current: 300 uA

Operating time: 7 minutes + / - 20 seconds

Warranty: 12 months

Features


convenient removable bracelets electrodes




powered by battery or adapter



suitable indicators and buttons




www.noraresearch.com
nora.research@mail.bg



Nora Research

Electronic device for bio-resonance therapy
(Hulda Clark zapper)



Tangra mini



Bio-resonance therapy - alternative

Purpose:

Electronic device "Tangra mini" is a new generation instrument is designed for active bio-resonance therapy (ABT) on the human body based on the generation of electromagnetic signals with certain frequencies and form by which the influence of harmful pathogens (viruses and bacteria) and parasites in human body. Typical for this type of therapy that is totally harmless, and that helps in the case where no other method does not help. Daily treatment with device clears the body of harmful parasites that can cause hundreds of diseases, while increasing tone, removes fatigue and improves overall health.

Made exactly according to original scheme of Dr. Hulda Clark zapper published in her book "Treatment of all diseases",

in this sense can be called "Hulda Clark Zapper". A scientist who laid the foundations of modern bio-resonance therapy through its practice and great books ("Treatment of all diseases" etc.). Apparatus is characterized by its comfort, quality and longevity, making it suitable for use in everyday life, medical and therapeutic centers. There are comfortable and healthy bracelets electrodes that make carrying comfortable therapy.

Have been used modern technologies such apparatus is designed and manufactured entirely in the European Union. Available is a detailed instruction manual in English. There is a guarantee and warranty support. More information can be found on our website:

www.noraresearch.com

Technical characteristics:

Power supply: 9V battery type 6F22

Power consumption: <0.015 A

Spectrum: a rectangular signal with a frequency of 32KHz (+ / - 10% when battery voltage 9V)

Amplitude of the output: max. 8V

Maximum output current: 0.006 A

Power consumption: max. 0.15W

Battery life: 200 cycles on average 7 minutes (depending on model of the battery)


Time of operation: manual mode

Mass Product Packaging: not more than 600 g

Warranty: 12 months

Features

convenient removable bracelets electrodes



Zapper in the original scheme of Dr. Hulda Clark



Comfortable and healthy jacks, compact design



www.noraresearch.com
nora.research@mail.bg



Nora Research

Electronic device to obtain a silver water
(colloidal silver generator)



Ponor



Silver ions - a source of health and life

Purpose:

Electronic device for silver ionization water "Ponor" is a new generation of equipment designed for the enrichment (dressing) of water (spring, tap, etc.) silver ions (Ag), thereby producing silvery water (colloidal silver). With his help can be continuously enriched random amounts of water. It is possible to obtain water at different concentrations, as well as for everyday use, preserving food and beverages and for higher concentrations of silver treatment and prolonged storage of water. The effect of silver ions is on more than 650 types of bacteria (as opposed to: the effects of any antibiotic is over no more than 5 to 10 types of bacteria). It is known that the antimicrobial effect of silver is higher than that of penicillin. Meanwhile, silver ions spare cells beneficial bacteria, absent adverse effects on the body

- something typical of antibiotics. Healing properties of colloidal silver are known from ancient times. Silver water can be used for everyday use for cosmetic and medicinal uses.

Device "Ponor" is characterized by its comfort, quality and long life, which makes it suitable for use in everyday life, medical and cosmetic centers. Silver electrode is made most pure silver and is intended for multi use. The use of the device is not bound by the use of a particular vessel. Modern microelectronics is used as the apparatus is designed and manufactured entirely in the European Union. Available is a detailed instruction manual in English. There is a guarantee and warranty support. For more information:

www.noraresearch.com

Technical characteristics:

Power: 220V/50Hz AC +/- 10% (optional 110V)

Power consumption: <0.5W

Amplitude of the output: 12V

Concentration of silver ions in accordance with the volume: 0.3 - 0.4 mg / l

Output current: reverse - 0.3Hz

Electrode resource: > 20 000 l (Resource of the electrode is approximate and applies only to wear the silver part, mechanical damage may occur much earlier);

Volume of ionized water: 1, 2 or 3 liters

Silver electrode: pure silver .999

Warranty: 12 months

Features


applicable to all vessels




electrode of pure silver with a long life



suitable indicators and buttons



www.noraresearch.com
nora.research@mail.bg



Nora Research

Silver ionizer for enrichment with silver ions and the decontamination of large volumes of water



Ponor pro



Silver ions - a source of health and life

Purpose:

Electronic device (apparatus) for silver ionization water "Ponor PRO" is a new generation of professional equipment designed for the enrichment (dressing) of water (spring, tap, etc.) silver ions (Ag), thereby producing Silver water (colloidal silver). Was designed specifically for high volume tanks (10 to 1000 liters), and wells and boreholes for drinking water daily consumption of 10 to 1000 liters. In use water is saturated with silver ions, while both purified from harmful microorganisms. Saving and drinking qualities for a long time, and acquires properties of a natural antibiotic. The effect of silver ions is on more than 650 types of bacteria (for comparison: the effect of any antibiotic is on no more than 5 to 10 types of bacteria). It is known that the antimicrobial effect of silver is higher than that of penicillin. At the same time the silver ions friendly bacteria cells, such as absent adverse effects on the body - something characteristic of

antibiotics. The healing properties of colloidal silver have been known since ancient times. The concentration of silver ions has guaranteed optimum values for drinking water. Apparatus is characterized by comfort, long life, stable parameters and automated operation combined with light control and functionality of the electrode. There are 3 programs work by the amount of water that can optionally be reprogrammed for different volumes (up to 1000 liters) on request. It has a silver electrode with increased weight made of silver with high quality (sample 0.999), with its cable length can reach 20 meters. Used is modern technology and precision microelectronics. Apparatus is designed and manufactured entirely in the European Union. Available detailed instructions in English. There is a warranty and post warranty maintenance. For more information:

www.noraresearch.com

Technical characteristics:

Power supply: 220V/50Hz AC +/- 10%
 Power consumption: <1W
 Amplitude of the output: 12V
 Output current: reversible with frequency 0.2 Hz
 Concentration of silver ions in accordance with the volumes: 0.02 to 0.05 mg/l
 Cable length: 5 m standard (option: up to 20m.)
 Electrode resource:> 30 000 l (electrode resource is approximate and only wear silver part, mechanical damage can occur much earlier);
 Volume of ionized water (standard):
 Program 1: from 10 to 50 liters
 Program 2: from 50 to 150 liters
 Program 3: from 150 to 300 liters
 (There are optional individual programs with other volumes with a maximum capacity of up to 1000 litres.)
 Electrode (anode and silver cathode) mass of silver: 10 g +/- 10%
 Mass Product Packaging: not more than 1.0 kg (up to 3.0 kg with optional extension cable).
 Warranty: 12 months

Features

Silver water in large quantities every day



Massive silver electrode with a long and robust cable



suitable indicators and buttons



www.noraresearch.com
nora.research@mail.bg



www.noraresearch.com
nora.research@mail.bg