



Nora Research

Professional electronic device for bio-resonance therapy
(professional Zapper)



Tangra II



Bio-resonance therapy - alternative

Purpose:

Professional electronic device "Tangra 2" is a new generation instrument is designed for active bio-resonance therapy (ABT) on the human body based on the generation of electromagnetic signals with certain frequencies and form by which the influence of harmful pathogens (viruses and bacteria) and parasites in human body. Typical for this type of therapy that is totally harmless, and that helps in the case where no other method does not help. Daily treatment with device clears the body of harmful parasites that can cause hundreds of diseases, while increasing tone, removes fatigue and improves overall health.

It is based entirely on the writings of the famous Dr. Hulda Clark, which laid the foundations of modern bioresonance therapy practice and through your great books ("Treatment of all diseases, etc..). There are 4 work therapy programs, which allows to have an influence on a wide range of harmful

microorganisms. One of the programs is that of electronic Zapper Dr. Hulda Clark. Available is a healing program for high-frequency therapy, which acts as an immune stimulant. Possible change in the power of the device, thereby achieving a better setting to the disease and given the sensitivity of human organism.

Device is characterized by its comfort, quality and long life, which makes it suitable for use in everyday life, medical and therapeutic centers. There are healthy and convenient removable bracelets electrodes. Power is carried through safe transformer adapter. Contemporary design, comfortable buttons, light and sound indicators make use of the device comfortably. Modern microelectronics is used as the unit. Is designed and manufactured entirely in the European Union. Available is a detailed instruction manual in English. There is a guarantee and warranty support. More information can be found on our website: www.noraresearch.com

Technical characteristics:

- Spectrum:

Program I: 30 000 Hz - rectangular shape with a wide harmonic range

Program II: 1 000 Hz - rectangular shape with a wide harmonic range

Program III: 5 000 Hz - rectangular shape with a wide harmonic range

Program Ultra: timing and frequency modulated signal with high frequency (HFT therapy) - a rectangular shape with a wide harmonic range

- **Output signal:** rectangular shape, 100% positive offset

- **Amplitude of the output:** max. 13V

- **Maximum output current:**

Mode "Low" 80 micro amps

Mode "High" 200 micro amps

- **Automatic timer** 7 minutes +/- 10%

- **Power of the device:** battery 9 volts/type 6F22 and optional version + through safety transformer power adapter from the mains 220/230V / 50Hz (optional 110V);

- **Power consumption:** <1W

- **Temperature range for normal operation:** from +10 to +35 degrees Celsius;


- **Operation in humidity:** <80%;

- **Weight:** not more than 2 kg;

- **Warranty:** 12 months (without transportation costs);

Features

convenient removable bracelets electrodes



4 programs, 2 power levels
suitable indicators and buttons



modern and comfortable design



www.noraresearch.com
nora.research@mail.bg